

Pluggo Advia Centaur 5 Position Automated Decapper - Technical Description

Pluggo is designed for the safe removal of original caps from most standard vacuum collection tubes. Two models are available adapted to handle 5ml or 7ml tubes, 75 x 13 mm or 100 x 13 mm size respectively, carried in Advia Centaur 5 position rack.

PLUGGO



A detachable chute guides the caps into a disposal area after decapping.

Technical data

Capacity:	Up to 5 samples in each rack
Sample tubes:	Standard vacuum blood tubes
Speed:	2 seconds per rack
Size:	29 x 27 x 21cm (Depth x Width x Height)
Weight:	8,0 kg
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA