


Pluggo Galileo – Technical description

Pluggo is designed for safe removal of original caps from BD Vacutainer™ vacuum collection tubes. Pluggo handles 5ml or 7ml Vacutainer™ vacuum tubes with 13 x 75 mm and 13 x 100 mm size tubes carried in Galileo 16-pos racks.

PLUGGO



The caps are guided to an outlet for collection of biohazard disposal with a detachable chute.

Technical data

Capacity:	Up to 16 samples in each rack
Sample tubes:	BD Vacutainer™
Speed:	Appox. One second per tube
Size:	435x500x210 mm (Depth x Width x Height)
Weight:	10,0 kg
Voltage:	230VAC50Hz/110VAC 60Hz
Power:	60VA