


## Pluggo RHs Hamilton 24 Position Racks Automated Decapper - Technical Description

Pluggo is designed for the safe removal of original caps from most standard vacuum collection tubes. Pluggo can handle 5ml, 7ml and 10ml vacuum tubes with 75 x 13 mm, 100 x 13 mm, 100 x 16 mm tubes carried in the Hamilton 24 position racks.

### PLUGGO



A detachable chute guides the caps into a disposal area beneath the instrument after decapping.

#### Technical data

Capacity:	Up to 24 samples in each rack
Sample tubes:	Standard vacuum blood tubes
Speed:	30 seconds per rack
Size:	42 x 70 x 21cm (Depth x Width x Height)
Weight:	9,0 kg
Voltage:	230VAC 50Hz / 110VAC 60Hz
Power:	60VA